



PT-MT3
Electro-Magnetic Test Tablet
Installation Qualification (IQ)
Version 5.2

Pharma Test Apparatebau AG
Siemensstrasse 5
D-63512 Hainburg (Germany)

Phone: (+49) 0 6182/9532-600
Fax: (+49) 0 6182/9532-80
E-Mail: info@pharma-test.de
Web: www.pharma-test.com



DIIN EN ISO 9001-Certificate No FS
529019/0388D

Table of Contents

Table of Contents	2
Document History.....	4
Introduction	5
Objective.....	5
Equipment.....	5
Instructions for Documentation Completion	6
Correcting Entries	6
Correcting Entries, sections and parts which are now required or available.....	6
Conditions Requiring Re-Qualification or re-calibration	6
Installation Qualification Program	7
Section 1.0 Instrument Identification.....	7
Section 2.0 Information User and Installation Place	7
Section 3.0 Required instrument, substances, and standards used for the Functional Qualification (OQ).....	7
Section 4.0 Documentation and Installation Procedure	7
Section 5.0 Instrument identification, supplied components, optional accessories.....	7
Section 6.0 Results and Comments.....	7
Section 1.0 Instrument Identification.....	8
Instrument Details.....	8
Section 1.1 Optional Accessories	8
Section 1.2 Required Documentation and Information.....	9
Ordering Information.....	9
Equipment Documentation	9
Section 1.3 Instrument Installation Site Requirements.....	9
Electrical Conditions at Installation Site	9
Section 1.4 Other Services Required.....	10
Section 1.5 Operating Environment.....	10
Section 1.6 Clearance	10
Section 2.0 Personnel Identification	11
Section 2.1 Place of Installation	12
Section 3.0 Required instrument, substances, and standards used for the Functional Qualification (OQ).....	13
Section 3.1 Certificates for the Standards, Instruments and Substances.....	13

Pharma Test Apparatebau AG

Installation Qualification (IQ)

Section 4.0	Documentation and Installation Procedure	14
Section 4.1	Unpacking.....	14
Section 5.0	Instrument identification, supplied components, optionals.....	15
Section 5.1	Mains Cord	15
Section 5.2	Complete the Instrument/System	15
Section 5.4	Connect Mains.....	16
Section 5.5	Connect the Measuring Head.....	16
Section 5.6	Turn on the Instrument	16
Section 5.7	Verify Firmware Version	16
Section 5.8	Turn off the Instrument	16
Section 5.9	Qualify the stop watch.....	17
Section 5.10	Qualify the Reference Weight PTB-CAL15	17
Section 5.11	Qualify the Reference Weight PTB-CAL30	17
Section 5.12	Qualify the Reference Weight Rise 2x10kg.....	17
Section 6.0	Result and Comments.....	18

Document History

Revision	Valid from	Author	Change	Remark
5.0	12.07.2012	PTAG	N	New document
5.1	13.08.2015	PTAG	R	PTB420 incl.
5.2	15.03.2017	PTAG	R	Design, PTB-M and weights incl.

Index Information - Change:

N = New Document

C = Correction

R = Revision

Introduction

Objective

Installation Qualification (IQ) is used to verify that the Pharma Test Instrument/System (or components where applicable) has been installed according to the manufacturer's requirements. The results are documented and explained as passed (OK) or failed (NOK), compared to the requirements laid down in this form

Equipment

The Pharma Test PT-MT3 electro-magnetic tablet is composed of:

- The PT-MT3 instrument
- Standard Supply Scope
- Optional Equipment and Accessories

The PT-MT3 represents a system for adjusting, calibrating, and qualifying the hardness range of Pharma Test and other brand tablet hardness testers. The **PT-MT3** instrument has an electronically controlled magnetically held "finger" made of stainless steel and a load cell attached to it. The finger is held in position by an electromagnet which is always powered with maximum current. The finger is equipped with a electronical load cell which can be calibrated to allow a maximum force of 500N. The strike arm or jaw of the hardness tester makes contact with the load cell at the finger and tries to "break" the finger away from the electromagnet holding it in place. Depending on the force setting, the force required to make this break can be set stepplessly to maximum 500 N. The load cell has to be calibrated against traceable standards to produce a useful series of equivalent tablet break points. The PT-MT3 can be supplied fully certified and includes a built-in thermo printer

Instructions for Documentation Completion

All performers and reviewers must complete qualification forms using the following guidelines:

1. Complete all items on a form in full.
2. Document any deviation from defined protocols and expected results. Owner approval of protocol deviations must be documented before final approval signatures can be obtained.
3. Write additional comments on an addendum sheet when there is not enough space on a form to accommodate all comments. Follow these three steps when adding an addendum sheet:
 4. Use your initials to confirm the entry
 5. Enter the date of the additions (addendum)
 6. Number the addendum pages numerically
7. Insert the addendum sheet behind the original page
8. Make all entries in permanent ink.

Correcting Entries

If you need to make corrections on a form, use the procedures described below:

Correcting Entries, sections and parts which are now required or available

It is possible that certain information or requirements are not available or necessary for the instrument to be qualified. This information may be a full section, a part of it or procedure. Mark this element clearly, so that it is understood that it is not necessary in this case.

To correct a long entry or information block on a form follow this procedure:

1. Draw a diagonal line, through the wrong, invalid or incorrect information
2. Write the correction on a separate addendum sheet
3. Give a brief explanation of the change.
4. Sign it using your initial
5. Enter the date of the change

NOTE: All original entries must remain legible after any corrections have been made.

Conditions Requiring Re-Qualification or re-calibration

CAUTION: The following conditions require re-qualification or calibration:

- When a system modification has been completed which affects the installation qualification
- When this system is being removed from where it was originally installed

Installation Qualification Program

This document is divided in sections.

Section 1.0 Instrument Identification

In this section the instrument including details and serial number(s), as well as the required documentation are identified.

Section 2.0 Information User and Installation Place

This section the user of the equipment and installation place are defined

Section 3.0 Required instrument, substances, and standards used for the Functional Qualification (OQ)

This section includes information about the equipment, tools and substances required to perform the Functional Qualification (OQ), as well as information of the place of installation.

Section 4.0 Documentation and Installation Procedure

This section contains the installation procedure, test protocols and test results in a pass/fail format for each test.

Section 5.0 Instrument identification, supplied components, optional accessories

This section is used to control the completeness of the supply scope and optional accessories which are needed for the later use of the instrument/system

Section 6.0 Results and Comments

This section is used to document the result of the installation and for comments regarding the installation procedure.

Section 1.0 Instrument Identification

Check if the instrument/system according to the order confirmation has been received and enter the serial number and the type name. The serial number is printed on the type plate on the back of the instrument:

Part-No.	Instrument Description	Type	Rec.	Miss.	NA	Serial No.
29-18050	Electro-Magnetic Tablet	PT-MT3				

Instrument Details

Voltage

100/115 V ☐ or 230/240 V ☐ 50/60 cycl.

Firmware Version

Asset No. or Lab ID No.

Section 1.1 Optional Accessories

Part-No.	Instrument Description	Rec.	Miss.	NA	Serial No.
29-18001	Support for: PTB 111E/ER, PTB 111-500E/ER, PTB 311E, 411E, 511E				
28-18001	Support for: PTB 301, 302, 502				
29-18002	Support for: PTB-M				
29-03150	Support for: PTB 420, 420-500				
27-18001	Support for PTBA 211E, 211-500E				
80-18001	Support for WHT M, WHT ME				
29-18005	Support for: Sotax/Pharmatron 4M,8M				

Performed by: _____

Signature

Date: _____

DD/MM/YYYY

Section 1.2 Required Documentation and Information

Ordering Information

Customer Order Number

Supplier Delivery-Note or Invoice Number

Equipment Documentation

Part-No.	Description	Rec.	Miss.	NA
B-29-18050	Operating Manual			
IQ-29-18050	Installation Qualification			
OQ-29-18050	Operation Qualification			
NA	Instrument Log Book			
QC-29-28050	QC/DQ Test Report			
NA	Delivery Note			

Section 1.3 Instrument Installation Site Requirements

NOTE: This instrument requires electrical installation with a quality, noise free ground (earth). In this case, as in all electrical installations where liquid handling is an integral part of the instrument function, all mains (line) supplies should be protected with an RCD (Residual Current Detector) which typically trips out when current leakage to ground (earth) of approx. 30mA is detected.

Electrical Conditions at Installation Site

Voltage: 100/115V ☐ 230/240V ☐ Frequency: 60cycl. ☐ 50cycl. ☐

Performed by: _____

Signature

Date: _____

DD/MM/YYYY

Section 1.4 Other Services Required

The Instrument has to be installed on a level surface. The supporting bench must be able to support the weight of the Instrument, as well as other equipment installed on the same bench continuously.

NOTE: Pharma Test is not responsible for adequate site planning and installation. We are however, happy to advise if our help is requested.

Section 1.5 Operating Environment

Temperature:	15-40°C	Rel. Humidity:	20-80%
Physical Site:	Clean, dry and level bench	Weight Support:	Approx. 60 kg/m ²
Required Bench Space (Length/Depth/ Height), without hardness test instrument:	20 x 30 x 20 cm		

Section 1.6 Clearance

A minimum of 10cm of the rear and sides of the instrument is required.

Performed by:	_____	Date:	_____
	Signature		DD/MM/YYYY

Section 2.0 Personnel Identification

Installation Engineer (1):

Name (print)

Initials

Signature

Date (DD/MM/YYYY)

Installation Engineer (2):
(optional)

Name (print)

Initials

Signature

Date (DD/MM/YYYY)

Installation Engineer (2):
(optional)

Name (print)

Initials

Signature

Date (DD/MM/YYYY)

Released by:

Name (print)

Initials

Signature

Date (DD/MM/YYYY)

Performed by:

Signature

Date:

DD/MM/YYYY

Section 2.1 Place of Installation

Company Name:

Address:

Department/Bldg.:

Location/Bldg.:

Contact:

Telephone:

Fax:

E-Mail:

Performed by:

Signature

Date:

DD/MM/YYYY

Section 3.0 Required instrument, substances, and standards used for the Functional Qualification (OQ)

Part-No.	Description	Serial No.	Calibrated Until	OK	NA
10-61000	Stop Watch				
28-00290	15 kg Weight				
28-00300	30 kg Weight				
28-00350	PTB-CAL weight raise set 2x10kg for PTB CAL30 - total weight 50kg				

Section 3.1 Certificates for the Standards, Instruments and Substances

PartNo.	Description	Calibration Date	Rec.	Miss.	NA
101-3200	DKD certificate for stop watch				
281-2890	DKD certificate for PTB-CAL15				
281-2830	DKD certificate for PTB-CAL30				
281-2840	DKD certificate for 2x10kg weight raise				

Performed by: _____

Signature

Date: _____

DD/MM/YYYY

Section 4.0 Documentation and Installation Procedure

This Section describes the performance and documentation of the Installation Qualification of the instrument/system. Complete each section as described.

Full information about the use and operation of the instrument/system is available from the supplied Operating Manual

Section 4.1 Unpacking

Unpack the instrument. Make sure that there are no more parts inside the packing material or box. Make sure that the instrument has not been damaged during transportation

OK	NOK	NA

Performed by: _____ Date: _____
Signature DD/MM/YYYY

Section 5.0 Instrument identification, supplied components, optional

All parts required for the standard operation of the instrument are listed at the supplied delivery note (invoice). Each part is labeled with its part-no. Please check the completeness of the delivery scope using the delivery-note (or invoice) attached.

Description	OK	NOK
All items as per attached delivery note have been fully supplied		

In case any item is missing, please be so kind and inform Pharma Test AG immediately. Please mention your Purchase Order No., the instrument serial number, the part number which is missing, as well as the number of the delivery note. If the missed part is essential for the operation of the instrument, you can only proceed with the Installation Qualification after the part has been received.

Section 5.1 Mains Cord

These parts are included in the standard supply scope of the Instrument

Part-No.	No.	Description	Rec.	Miss.	NA
34-08500	1	IEC/EUR mains cord			
34-08511	1	US mains cord			
34-08512	1	GB mains cord			
34-08513	1	ARG/AUS/NZ mains cord			

Section 5.2 Complete the Instrument/System

Install and complete the instrument with its components as per supplied operating manual

OK	NOK	NA

Performed by: _____

Signature

Date: _____

DD/MM/YYYY

Section 5.4 Connect Mains

Keep main power switch of the instrument at off position. Connect and verify connection of the power cable into the ground fault protection AC outlet. Check that the mains cord and socket do not show any visible damage.

OK	NOK	NA

Section 5.5 Connect the Measuring Head

Use the supplied cable to connect the measuring head to the main instrument. Plug the cable into the corresponding port on the back side of the main PT-MT3 instrument.

pass	Fail	N/A

Section 5.6 Turn on the Instrument

Turn on the instrument following the description of the manual. Use the mains power switch at the rear side of the instrument. All displays will light up.

OK	NOK	NA

Section 5.7 Verify Firmware Version

After the instrument is turned on some instrument show the installed Firmware Version Number. Check that the displayed version number is the same as noted in the QC/DQ document for this instrument.

OK	NOK	NA

Section 5.8 Turn off the Instrument

Turn off the instrument using the power switch on the rear left side of the instrument.

OK	NOK	NA

Performed by: _____

Signature

Date: _____

DD/MM/YYYY

Section 5.9 Qualify the stop watch

Each instrument, which is used for the functional qualification (OQ) of this instrument has been tested and qualified at factory. A certificate issued from an authorized laboratory is only supplied when it has been ordered

The stop watch is used to qualify the time counter

OK	NOK	NA

Section 5.10 Qualify the Reference Weight PTB-CAL15

Each reference weight, which is used for the functional qualification (OQ) and calibration of this instrument has been tested and qualified at factory. A certificate issued from an authorized laboratory is only supplied when it has been ordered

The PTB-CAL15 is used to qualify and calibrate the digital load cell of the instrument: 15kg ± 15g (not calibrated) - 15kg ± 0.5g (calibrated incl. certificate)

OK	NOK	NA

Section 5.11 Qualify the Reference Weight PTB-CAL30

Each reference weight, which is used for the functional qualification (OQ) and calibration of this instrument has been tested and qualified at factory. A certificate issued from an authorized laboratory is only supplied when it has been ordered

The PTB-CAL30 is used to qualify and calibrate the digital load cell of the instrument: 30kg ± 30g (not calibrated) - 30kg ± 1g (calibrated incl. certificate)

OK	NOK	NA

Section 5.12 Qualify the Reference Weight Rise 2x10kg

Each reference weight, which is used for the functional qualification (OQ) and calibration of this instrument has been tested and qualified at factory. A certificate issued from an authorized laboratory is only supplied when it has been ordered

The Weight Raise 2x10kg is used to qualify and calibrate the digital load cell of the instrument: 50kg ± 50g (not calibrated) - 50kg ± 2.5g (calibrated incl. certificate)

OK	NOK	NA

Performed by: _____ Date: _____
Signature DD/MM/YYYY

**Pharma Test Apparatebau AG
Installation Qualification Testing
Certificate**

Pharma Test Apparatebau AG
Siemensstraße 5
D-63512 Hainburg/Germany

Telephone +49 6182 9532-600
Telefax +49 6182 9532-80
+49 6182 9532-650
info@pharma-test.de
www.pharma-test.com

Section 6.0 Result and Comments

Type		Serial Number	
Mains Voltage		Firmware Version	

The instrument has passed the Installation Qualification (IQ) procedure.

Yes ☐

No ☐

Check yes if all tests have passed. In case one or more tests failed check no and document the reason for the failure on this report. In this case the installation qualification has to be repeated once the reason for failure has been eliminated.

Comments

This completes the Installation Qualification of the Instrument.

Performed by: _____ Date: _____
Signature DD/MM/YYYY

Approved by: _____ Date: _____
Signature DD/MM/YYYY

Released by: _____ Date: _____
Signature DD/MM/YYYY